## ▲ **ALTEC**<sup>™</sup>

Email : sales@equipmentsexporters.com Phone: +91-9311469084

Product Name : Code :   Gel Developing Tray ALABS-A89-049	
Description :	
Gel Developing Tray Used for staining, fixation, destaining and handling of electrophoresis gels and membrances. It is made of ABS Plastic injection moulded with acrylic cover. It is tapered from one end to allow easy liquid drainage with minimum gel handling.	
Size:	
<b>L</b> = 142 mm	
<b>W</b> = 120 mm	

**H** = 58 mm

## **Equipments Exporters**

Website: www.equipmentsexporters.com, Email: sales@equipmentsexporters.com Address: 75, Lajpat Nagar-IV, New Delhi-110024 Phone: +91-9311469084

