

**Product Name :**  
Steam Distillation Set Up**Product Code :**  
ALABS-A104-073**Description :****Steam Distillation Set Up**

The equipment consists of jacketed pressure vessel for distillation. Steam from a steam generator can either be allowed to enter 'in the outer jacket or in the vessel. The vapors from the top of the vessel are condensed in the shell and tube type condenser by circulating cooling water supplied by laboratory overhead tank.

The condensate is collected in a separating chamber. The set-up is fitted with steam trap control valves pressure gauge, temperature sensors and other instrumentations required

**Experiment :**

- To study the Packed Distillation Column

**Feature :**

- Stainless Steel tanks and wetted parts
- Superb Painted structure
- Simple to operate & maintain
- Compact & stand alone set up

**Utilities Required :**

- Water Supply 2 LPM at 5 m head
- Drain
- Electricity Supply: 1Phase, 220 V AC, 6.5 kW
- Required Chemicals
- Refractrometer for analysis



## **Equipments Exporters**

**Website:** [www.equipmentsexporters.com](http://www.equipmentsexporters.com), **Email:** [sales@equipmentsexporters.com](mailto:sales@equipmentsexporters.com)

**Address:** 75, Lajpat Nagar-IV, New Delhi-110024 **Phone:** +91-9311469084