# ▲ ALTEC<sup>™</sup>

Email : sales@equipmentsexporters.com Phone: +91-9311469084

	Product Code :
Steam Distillation Set Up	ALABS-A153-01



#### **Description :**

Steam Distillation Set Up

The equipment consists of jacketed

pressure vessel for distillation. Steam from a steam generator can either be allowed to enter 'in the outer jacket or in the vessel. The vapors form the top of the vessel are condensed in the shell and tube type condenser by circulating cooling water supplied by laboratory overhead tank.

The condensate is collected in a separating chamber. The set-up is fitted with steam trap control valves pressure gauge, temperature sensors and other instrumentations required

### Experiment :

• To study the Packed Distillation Column

Feature :

- Stainless Steel tanks and wetted parts
- Superb Painted structure
- Simple to operate & maintain
- Compact & stand alone set up

#### **Utilities Required :**

- Water Supply 2 LPM at 5 m head
- Drain
- Electricity Supply: 1Phase, 220 V AC, 6.5 kW
- Required Chemicals
- Refractrometer for analysis

# ▲ ALTEC<sup>™</sup>

### **Equipments Exporters**

Website: www.equipmentsexporters.com, Email: sales@equipmentsexporters.com Address: 75, Lajpat Nagar-IV, New Delhi-110024 Phone: +91-9311469084