

Product Name :
Clinostat Electrical**Product Code :**
ALTECEDULABC1066**Description :****Clinostat Electrical**

For demonstration of heliotropism and geotropism. Germinating seedlings are pinned to the cork table which is rotated by an electric motor at 4 revolutions per hour. The motor is spindle mounted and may be clamped by the thumbscrew so that the plane of the table can be set at any angle between horizontal and vertical. An Acrylic cover is provided for the cork table to shield the seedling from draughts. Diameter of cork disc 100 mm.